**GenePharma**

# **DxGene™ Tissue and Cell Total RNA Extraction Kit 1**

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---

## Contents

<b>I. Overview</b>	
<b>Introduction and Principle</b> .....	11
<b>Kit Content and Specification</b> .....	11
<b>Equipment and Reagents to Be Supplied by User</b> .....	12
<b>Related Products</b> .....	13
<b>II. Protocol</b>	
<b>Sample Type Requiement</b> .....	14
<b>Total RNA isolation from Tissue and Cell Line</b> .....	14
<b>III. Troubleshooting</b>	
<b>Troubleshooting</b> .....	15
<b>IV. Downstream Applications</b>	
<b>Applications</b> .....	18
<b>Related Products and Service</b> .....	18
<b>V. Appendix</b>	
<b>Storage of RNA Samples</b> .....	19
<b>Quantitation of RNA</b> .....	19
<b>Judgement of the Qaulity of Total RNA</b> .....	20
<b>Technical Service</b> .....	22

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产品目录号	产品名	规格	储存条件
QPG-060	DxGene™ Tissue and Cell line Total RNA Extraction Kit	20 Samples	室温
QPG-061	DxGene™ Tissue and Cell line Total RNA Extraction Kit	50 Samples	室温

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## I. Overview

### Introduction and Principle

DxGene™ Tissue and Cell line Total RNA Extraction Kit is developed and optimized for total RNA of high quality extraction from cell lines and tissue samples including lipid tissue samples. It can be used with DxGene™ Total RNA Quantitation Kit, DxGene™ Her2 Test Kit, DxGene™ EGFR Test Kit, DxGene™ Beta-actin Test Kit and DxGene™ Her2 Duplex Test Kit, providing total RNA of high quality for detecting and quantitation of Her2 and EGFR gene accurately in cell and tissue samples

DxGene™ Tissue and Cell line Total RNA Extraction Kit integrates traditional phenol/guanidine-based lysis and silica-gel-membrane purification of total RNA. The kit also provides Ezol™ lysis reagent, designed to facilitate lysis of samples and inhibit RNase. After RNA binding to the silica membrane, phenol and other contaminations are efficiently washed away by the WB buffer in the kit. Then high quality of total RNA is eluted in 50µl or more RNase-free water.

### Kit Specification

Cat. No.	Item	Specificaion	Storage
QPG-060	DxGene™ Tissue and Cell line Total RNA Extraction Kit	20 Samples	RT
QPG-061	DxGene™ Tissue and Cell line Total RNA Extraction Kit	50 Samples	RT

---

**Kit Contents**

Contents	Specification	
	20 Samples	50 Samples
Mini-spin Column	20	50
1.5ml RNase-free tube	20	50
2.0ml Collection tube	20	50
Ezol™ Reagent	24ml	60ml
WB Buffer (Concentrate)	11ml	27.5ml
RNase-free Eution Buffer	2ml	5ml
User Manual	1	1

---

**Equipment and Reagents to Be Supplied by User**

---

When working with chemicals, always wear a lab coat, disposable gloves and protective goggles.

**Reagent**

Chloroform  
 95%~100% Ethanol  
 RNase-free water/DEPC water

**Equipment**

Sterile, RNase-free tips  
 Disposable gloves  
 Equipment and tubes for tissue disruption and homogenization.  
 RNase-free 1.5ml or 2ml microcentrifuge tubes  
 Microcentrifuges (with a rotor for 2ml tubes)  
 for centrifugation at 4°C and room temperature.

---

**Related Products**

DxGene™ Total RNA Quantitation Kit, DxGene™ Her2 Test Kit, DxGene™ EGFR Test Kit, DxGene™ Beta-actin Test Kit and DxGene™ Her2 Duplex Test Kit are also provided by GenePharma. You can contact us for more information.

<b>Cat. No.</b>	<b>Item</b>	<b>Specification</b>
QPG-070	DxGene™ Total RNA Quantitation Kit	20 rxns
QPG-080	DxGene™ Her2 Test Kit	100rxns
QPG-081	DxGene™ EGFR Test Kit	100rxns
QPG-082	DxGene™ Beta-actin Test Kit	100rxns
QPG-083	DxGene™ Her2 Duplex Test Kit	100rxns

---

## II. Protocol

---

### Sample type requirement

DxGene™ Tissue and Cell line Total RNA Extraction Kit is developed and optimized for total RNA isolation in fresh tissue samples or cell samples. As to RNA extraction in FFPE samples, please select DxGene™ FFPE Total RNA Extraction Kit.

---

### Protocol

1. To extract RNA from tissues, cut a piece of deep-frozen tissue in a chilled mortar, and weight by electronic balance. Transfer tissue immediately into a 1.5-ml RNase-free tube which already contains 1 ml of Ezol reagent. Homogenize tissue with an electronic homogenizer on ice for 15 s or until no visible tissue, and then move to step 3.

*Note: Tissue sample should be stored in liquid nitrogen or in deep freezer (-80°C). This protocol is suitable for isolation of RNA from less than 50 mg of tissue; otherwise, more Ezol reagent is needed to digest tissue sample.*

2. To extract RNA from cells, put cell culture dish or flask on ice, remove culture medium, and wash cells twice with ice-cold 1X PBS. Add 1 ml of Ezol reagent and scrape cells off the plate with cell lifter. Transfer cell lysates into 1.5-ml RNase-free tube, vortex thoroughly for 30 s, and then move to step 3.

*Note: PBS buffer is not provided with kit, and is available separately.*

*Note: For suspension cells, spin down cells first at 2,000 x g for 3 min, wash cells twice with 1X PBS, and then add 1 ml of Ezol reagent.*

*Note: This protocol is suitable for isolation of RNA from  $\leq 1 \times 10^7$  cells, or from a  $\leq 10$ -cm diameter dish; otherwise, more Ezol reagent is needed to homogenize cells.*

3. Incubate mixture at room temperature for 5 min, and add 0.2 ml of chloroform. Cap the tube securely, and shake it vigorously for 15 s.

*Note: Avoid vortexing as this may increase the DNA contamination to the RNA sample.*

4. Place the tube at room temperature for another 2-3 min, and centrifuge at 12,000 x g for 15 min at 4 °C.

*Note: Centrifugation separates solution into three phases: upper aqueous with RNA, middle layer with protein/DNA, bottom organic with DNA. The 4 °C is essential for phase separation.*

5. Transfer the upper aqueous phase to a new 1.5-ml sterile tube, and add 1 volume of 70% ethanol. Mix solution thoroughly by vortex for 5 s.

*Note: Discard the middle layer and bottom organic phase if no need for isolation of DNA.*

6. Pipet up to 700  $\mu$ l of the sample immediately, including any precipitate that may have formed, into a mini-spin column in a 2-ml collection tube. Close tube gently, and centrifuge at 8,000 x g for of 15 s at room temperature. Discard the flow-through.
  7. Repeat step 6 using the remaining sample, and discard the flow-through.
-

8. Add 700  $\mu$ l of Wash Buffer (WB) into a mini-spin column. Close the tube gently, and centrifuge at 8,000  $\times$   $g$  for 15 s. Discard the flow-through.

*Note: Ensure that three volume of 100% ethanol is added to WB before using for the first time.*

9. Repeat step 8 twice with 500  $\mu$ l wash buffer. After last wash, centrifuge column for 2 min at 8,000  $\times$   $g$  to dry silica-gel membrane.

*Note: It is important to completely remove washing solution. Any carryover may affect elution.*

10. Transfer column into a new RNase/DNase-free 1.5-ml tube, pipet 40  $\mu$ l RNase/DNase-free water directly onto the RNeasy silica-gel membrane. Close the tube gently, wait 1 to 5 min, and centrifuge for 1 min at 8,000  $\times$   $g$  to elute RNA.

11. Load elute back to the column, and repeat the elution step 10.

12. Keep eluted RNA sample on ice all the time and store at  $-70^{\circ}\text{C}$ .

### III. Troubleshooting

---

Kit is carefully assembled and contains high quality of reagents. This troubleshooting is to provide a general guide for solving the problems which are often related to tissue handling and operation. Many other problems may occur during extraction. You are always welcome to discuss any problem with our technique service.

---

#### Phases don't separate completely

##### Insufficient homogenate and mix

This is usually caused by insufficient mix after addition of chloroform, or by centrifugation at a high temperature. As described in step 3 and 4 of the protocol, you need to shake the tube vigorously for at least 15 s after addition of chloroform. Centrifuge the tube at 4°C or no higher than 8°C. A higher temperature will disrupt phase separation.

---

#### Clogged Column

##### Insufficient homogenate

Insufficient homogenate may result in large tissue piece that may clog the column. Please increase the g-force and time of centrifugation if necessary.

##### Too much starting material

The starting material can not be more than 100mg tissue. Please reduce the amount of starting material if necessary.

---

#### Degradation of RNA

##### Sample or material stored in bad condition

Total RNA in cells or tissues that in bad condition or dealed with some medication may lead to self-degradation. Please store samples in -80C or liquid nitrogen. It should be better if some commercial RNase inhibitor such as RNA Later.

##### Sample Collection

It is important to inactivate RNase as soon as possible once samples are gathered. Please shorten the time of collection as possible as you can to yield higher quality of total RNA.

**Sample manipulation before lysis** Keep tissue frozen before lysis. Please transfer tissue sample into lysis buffer contains RNase inhibitor before thawing and homogenize it in ice bath immediately.

**Centrifugation** Please centrifuge samples at 4C after adding chloroform.

**Degradation after electrophoresis** If loading RNA is more than 5µg, the band will show some smear. So please take about 700ng to 800ng total RNA for electrophoresis.

**Exogenous RNase contamination** Disposable gloves and respirator must be used during the whole process to prevent any RNase contamination. The mortar and muller should be dried and heated at 180°C over night before use to eliminate RNase. All the tubes and tips needed should be treated by DEPC water and then autoclave them before use.

RNA is degraded by ribonucleases (RNases) which are very stable and active enzymes. Although the materials and reagents provided by kit are free of RNases, RNases may be present in the working environment, such as on hands, laboratory bench, and equipment. Therefore, it is perceivable that RNA samples could be contaminated by RNases during or after isolation procedure. We suggest you take following precautions to avoid inadvertent RNase contamination. **1)** Use RNase-free and disposable supplies; **2)** Clean up the bench top or equipment with RNase inhibitor-containing solution; **3)** Always wear latex or vinyl gloves while handing reagents and RNA samples; **4)** Collect tissue samples as quickly as possible and avoid contact with any surface; **5)** Keep purified RNA samples on ice while preparing for downstream analysis; **6)** Store RNA sample at -70°C in aliquot to reduce freeze/thaw cycle; **7)** Clean electrophoresis tank with detergent solution (e.g., 0.5% SDS) and use DEPC-treated water to prevent RNA degradation in gel running; **8)** Always dilute RNA samples with RNase/DNase-free water provided with kit.

**Low Yield** It is normal to see the variation of RNA yield. However, an extreme low yield of RNA may indicate failure of extraction. Many factors can reduce the yield of RNA, such as poor quality of tissue samples, insufficient homogenization, and poor elution. It is important to select good quality tissue and handle it appropriately. Complete homogenization in Ezol reagent is also essential for extraction of RNA. In some

cases, a low yield of RNA is caused by poor elution, and an extended incubation with RNase-free water on column may enhance RNA recovery.

---

**Low  $A_{260}/A_{280}$  Value**

Usually, a low  $A_{260}/A_{280}$  value is caused by measurement of absorbance in water, or by contamination of proteins, phenol, or other organic chemicals. Under normal conditions, the isolated RNA should be free of phenol and organic chemicals. However, insufficient homogenization may cause co-purification of proteins and increase reading at 280 nm. We recommend following the described procedures to digest and homogenize tissue samples. We also recommend measuring the ratio of  $A_{260}/A_{280}$  in 10 mM Tris-HCl buffer.

## IV. Downstream Application

---

### Application

Isolated total RNA can be applied for many kinds of experiments such as real-time PCR, northern blotting or genechip that you need. The high purity and quantity RNA extracted may provide you more consistent result. Besides, this kit also retain small RNA including microRNA that you can also detect for research.

---

### Related products and Service

Also, for getting consistent results for gene detection, this kit is sold paired with DxGene™ Total RNA Quantitation kit, DxGene™ Her2 Testing Kit, DxGene™ EGFR Testing Kit and DxGene™ Duplex Testing Kit. It is convenient for researcher to deal with samples from extraction, quantitation to the detection. The technical support in Genepharma is always happy to answer any questions you may have either the information or the protocol in this manual and other related products. Please see the appendix for more contact information.

---

## V. Appendix

---

### Storage of RNA Samples

RNA is eluted in RNase/DNase-free water. There is no need to add any other reagents, such as RNase inhibitors or stabilizers. The eluted RNA samples can be stored at -20 or -70°C. Storage at -70°C is recommended for long-term usage. Under this condition, no obvious degradation of RNA is detectable after 3 years. To avoid possible degradation caused by frequent freeze/thaw, samples may be aliquoted before storage.

---

### Quantitation of RNA

RNA samples prepared by this kit are basically free of contamination by any organic chemicals, and therefore, are well suitable for quantitation by spectrophotometer method. The concentration of RNA can be determined by measuring the absorbance at 260 nm ( $A_{260}$ ) using quartz or UV-compatible cuvettes. An absorbance of 1 unit at 260 nm is equivalent to 40  $\mu\text{g}$  of RNA per ml. This relationship is valid only in water, so RNA samples should be diluted or measured in RNase/DNase-free water. To minimize variation, readings of sample absorbance should be between 0.1 and 1.0. Given that the isolated RNA is usually in small quantity, we recommend using ultra-micro cuvettes, which consume only one tenth of RNA sample required for measurement in regular cuvettes. UV compatible ultra-micro cuvettes are available from several biotech companies, such as BrandTech Scientific ([www.brandtech.com](http://www.brandtech.com)) and Santa Cruz Biotechnology ([www.scbt.com](http://www.scbt.com)). Before measuring samples, we also recommend calibrating spectrophotometer with RNA standard solution, which is available separately.

**Note: An example of determination of RNA concentration using regular cuvette:**

*Total volume of RNA sample, 50  $\mu\text{l}$*

*Dilution Factor, 1:100 (10  $\mu\text{l}$  of RNA sample + 990  $\mu\text{l}$  of sterile water)*

*$A_{260}$  Reading, 0.20 (measured in a 1-ml cuvette)*

*So, Concentration of RNA sample =  $40 \times A_{260} \times \text{dilution factor}$*   
*=  $40 \times 0.20 \times 100 = 800 \mu\text{g/ml}$*

*And, Total yield of RNA sample =  $800 \times 0.05 = 40 \mu\text{g}$*

---

***Determination of RNA concentration for the same sample using ultra-micro cuvette:***

*Total volume of RNA sample, 50  $\mu$ l*

*Dilution Factor, 1:100 (1  $\mu$ l of RNA sample + 99  $\mu$ l of sterile water)*

*A<sub>260</sub> Reading, 0.20 (measured in an ultra-micro cuvette)*

*So, Concentration of RNA sample = 40 x A<sub>260</sub> x dilution factor*  
*= 40 x 0.20 x 100 = 800  $\mu$ g/ml*

*And, Total yield of RNA sample = 800 x 0.05 = 40  $\mu$ g*

Customers may also use other methods, such as RiboGreen (Invitrogen) or NanoDrop equipment (Thermo Scientific), to determine RNA concentration for the samples prepared by this kit, and should obtain comparable results as spectrophotometer described above. However, the customers need to prepare a RNA standard curve in RiboGreen assay, or to calibrate NanoDrop machine by a RNA standard before measuring samples. Customers are also encouraged to follow manufactory's manual if other methods are used for quantitation of RNA samples.

---

**Judgement of the Quality of Total RNA****Integrity**

The integrity or the size of total RNA purified by this kit can be checked by denaturing agarose gel electrophoresis, or by a Bioanalyzer. Normally, there will be two sharp bands on gel following ethidium bromide staining, which represent 18S rRNA and 28S rRNA. If a sample shows good integrity of RNA, the ratio of the intensity of 28S and 18S will be close to 2:1. Bioanalyzer can provide more accurate analysis of RNA integrity.

***Note: Given that GeneMedDx™ real-time PCR kits are designed to amplify the mRNA fragments of shorter than 100 nucleotides, a compromised RNA integrity is expected to have marginal impact on the real-time PCR assay.***

**Purity**

This kit has integrated classic guanidine/phenol/chloroform extraction method and column purification process, and should produce high purity samples under the conditions described in the kit. Therefore, we do not recommend testing the purity. If, for some reasons, customers need to document sample purity, the ratio of the readings at 260 nm and 280 nm ( $A_{260}/A_{280}$ ) can be measured to provide an estimate of the purity of RNA.

Generally speaking, a pure RNA sample has an  $A_{260}/A_{280}$  ratio of 1.9-2.1. However, the  $A_{260}/A_{280}$  ratio is influenced considerably by pH. Since water is not buffered, the  $A_{260}/A_{280}$  ratio measured in water is often lower than that in buffer solution, and may vary greatly. As a consequence, the sample may be misinterpreted to have been contaminated with protein, phenol or other organic chemicals which usually have high absorption at 280 nm. For accuracy, we recommend measuring absorbance in 10 mM Tris-HCl (pH 7.5), and be sure to calibrate spectrophotometer with the same buffer.

*Note: RNA concentration should be determined in water, since the relationship between absorbance and concentration is based on an extinction coefficient calculated for RNA in water.*

*Note: This kit is optimized to prevent DNA contamination. Our tests indicate that there is no need to treat the samples with DNase. RNA samples prepared by this kit are fully compatible with GeneMedDx™ real-time PCR kits, which are designed to have the primers and probes specific for amplification of mRNA molecules, and also to contain 'No RT' control in the assay for monitoring possible interference by DNA.*

---

## Technical Service

**World Wide Web** Visit the GenePharma Web Resource using your World Wide Web browser. At the site, you can:

- Explore our catalog with full color graphics
- Get the scoop on our hot new products and special product offers
- Request catalog and product literature

The GenePharma URL is [www.genepharma.com](http://www.genepharma.com)

---

## Quality Control

The product has good QC system and is tested functionally by qRT-PCR using both MCF7 and MDA-MB-231 total RNA as template for every batch.

---

## Contact Us

For more information or technical assistance, please call, write, fax, or email. Additional international offices are listed on our web page ([www.genepharma.com](http://www.genepharma.com)).

### Corporate Headquarters:

Genepharma Corporation  
1011 Halley Road,  
Zhangjiang High-tech Park,  
ShangHai, CHINA  
Tel: 86-21-51320195  
Fax: 86-21-51320295  
E-mail: [service@genepharma.com](mailto:service@genepharma.com)